

Teklad Global 19% Protein Extruded Rodent Diet

Product Description- 2019 is a fixed formula, non-autoclavable diet manufactured with high quality ingredients and designed to support gestation, lactation, and growth of rodents. It is particularly useful for poorly breeding strains of rodents and genetically engineered stocks and strains. 2019 does not contain alfalfa or soybean meal, minimizing the occurrence of natural phytoestrogens, which can be beneficial for reproductive toxicology. Typical isoflavone concentrations (daidzein + genistein aglycone equivalents) range from non-detectable to 20 mg/kg. Exclusion of alfalfa reduces chlorophyll, improving optical imaging clarity. Absence of animal protein and fish meal minimizes the presence of nitrosamines. The extruded form ensures that despite its high fat content, it remains firm with minimal wastage. Also available irradiated (2919). For autoclavable diet, refer to 2019S (Sterilizable).

| Macronutrionts | , | |
|--------------------------------------|---------------|------------|
| Macronutrients | 0/ | 10.0 |
| Crude Protein | % | 19.0 |
| Fat (acid hydrolysis) a | % | 9.0 |
| Carbohydrate (available) b | % | 44.9 |
| Crude Fiber | % | 2.6 |
| Neutral Detergent Fiber ^c | % | 12.1 |
| Ash | % | 5.0 |
| Energy Density ^d | kcal/g (kJ/g) | 3.3 (13.8) |
| Calories from Protein | % | 23 |
| Calories from Fat | % | 22 |
| Calories from Carbohydrate | % | 55 |
| Minerals | | |
| Calcium | % | 0.9 |
| Phosphorus | % | 0.7 |
| Non-Phytate Phosphorus | % | 0.4 |
| Sodium | % | 0.1 |
| Potassium | % | 0.4 |
| Chloride | % | 0.4 |
| Magnesium | % | 0.2 |
| Zinc | mg/kg | 60 |
| Manganese | mg/kg | 80 |
| Copper | mg/kg | 15 |
| lodine | mg/kg | 6 |
| Iron | mg/kg | 200 |
| Selenium | mg/kg | 0.23 |
| Amino Acids | 3 3 | |
| Aspartic Acid | % | 1.1 |
| Glutamic Acid | % | 3.5 |
| Alanine | % | 1.2 |
| Glycine | % | 0.7 |
| Threonine | % | 0.6 |
| Proline | % | 1.8 |
| Serine | % | 0.9 |
| Leucine | % | 2.3 |
| Isoleucine | % | 0.8 |
| Valine | % | 0.9 |
| Phenylalanine | % | 1.1 |
| | | 0.5 |
| Tyrosine Methionine | | 0.5 |
| | | |
| Cystine | % | 0.3 |
| Lysine | % | 0.9 |
| Histidine | % | 0.4 |
| Arginine | % | 0.8 |
| Tryptophan | % | 0.2 |

Teklad Diets are designed and manufactured for research purposes only.



Ingredients (in descending order of inclusion)- Ground wheat, ground corn, corn gluten meal, wheat middlings, soybean oil, calcium carbonate, dicalcium phosphate, brewers dried yeast, L-lysine, iodized salt, magnesium oxide, choline chloride, DL-methionine, calcium propionate, L-tryptophan, vitamin E acetate, menadione sodium bisulfite complex (source of vitamin K activity), manganous oxide, ferrous sulfate, zinc oxide, niacin, calcium pantothenate, copper sulfate, pyridoxine hydrochloride, riboflavin, thiamin mononitrate, vitamin A acetate, calcium iodate, vitamin B₁₂ supplement, folic acid, biotin, vitamin D₃ supplement, cobalt carbonate.

| Standard | Product Form: | Extruded |
|--|---------------|----------|
| Vitamins | | |
| Vitamin A ^{e, f} | IU/g | 15.0 |
| Vitamin D ₃ e, g | IU/g | 1.5 |
| Vitamin E | IU/kg | 110 |
| Vitamin K ₃ (menadione) | mg/kg | 50 |
| Vitamin B ₁ (thiamin) | mg/kg | 17 |
| Vitamin B ₂ (riboflavin) | mg/kg | 15 |
| Niacin (nicotinic acid) | mg/kg | 75 |
| Vitamin B ₆ (pyridoxine) | mg/kg | 18 |
| Pantothenic Acid | mg/kg | 33 |
| Vitamin B ₁₂ (cyanocobalamin) | mg/kg | 0.08 |
| Biotin | mg/kg | 0.40 |
| Folate | mg/kg | 4 |
| Choline | mg/kg | 1200 |
| Fatty Acids | | |
| C16:0 Palmitic | % | 0.9 |
| C18:0 Stearic | % | 0.2 |
| C18:1ω9 Oleic | % | 1.7 |
| C18:2ω6 Linoleic | % | 3.9 |
| C18:3ω3 Linolenic | % | 0.4 |
| Total Saturated | % | 1.2 |
| Total Monounsaturated | % | 1.7 |
| Total Polyunsaturated | % | 4.4 |
| Other | | |
| Cholesterol | mg/kg | |

^a Ether extract is used to measure fat in pelleted diets, while an acid hydrolysis method is required to recover fat in extruded diets. Compared to ether extract, the fat value for acid hydrolysis will be approximately 1% point higher.

For nutrients not listed, insufficient data is available to quantify.

Nutrient data represent the best information available, calculated from published values and direct analytical testing of raw materials and finished product. Nutrient values may vary due to the natural variations in the ingredients, analysis, and effects of processing.

^b Carbohydrate (available) is calculated by subtracting neutral detergent fiber from total carbohydrates.

^c Neutral detergent fiber is an estimate of insoluble fiber, including cellulose, hemicellulose, and lignin. Crude fiber methodology underestimates total fiber.

d Energy density is a calculated estimate of metabolizable energy based on the Atwater factors assigning 4 kcal/g to protein, 9 kcal/g to fat, and 4 kcal/g to available carbohydrate.

^e Indicates added amount but does not account for contribution from other ingredients.

f 1 IU vitamin A = 0.3 µg retinol

g 1 IU vitamin D = 25 ng cholecalciferol