TD.94047 Vitamin Mix, AIN-93-VX

Formula	g/Kg
Niacin	3.0
Calcium Pantothenate	1.6
Pyridoxine HCl	0.7
Thiamin HCI	0.6
Riboflavin	0.6
Folic Acid	0.2
Biotin	0.02
Vitamin B ₁₂ (0.1% in mannitol)	2.5
Vitamin E, DL-alpha tocopheryl acetate (500 IU/g)	15.0
Vitamin A Palmitate (500,000 IU/g)	0.8
Vitamin D ₃ , cholecalciferol (500,000 IU/g)	0.2
Vitamin K ₁ , phylloquinone	0.075
Sucrose, fine ground	974.705

Reference: J. Nutr. 123:1939-1951, 1993. Supplies the recommended concentrations

of vitamins for AIN-93G and AIN-93M diets, when used at 10 g/kg of diet. The choline

source should be listed as a separate item in diets utilizing this vitamin mix.

Teklad Mixes are designed & manufactured for research purposes only.

Key Features

- Vitamin Mix
- Commonly Used

Key Planning Information

- Products are made fresh to order
- Lead time 2 weeks
- Store product at 4°C or lower
- Use within 6 months (applicable to most mixes)
- Box labeled with product name, manufacturing date, and lot number



Product Specific Information

- Powder
- Minimum order 500 gms

Options (Fees Will Apply)

- Rush order (pending availability)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- · Outside U.S.A or Canada ·
- askanutritionist@harlan.com

Speak With A Nutritionist

(800) 483-5523

Footnote

askanutritionist@harlan.com

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220 www.harlan.com



Place Your Order (U.S.A & Canada)

- Place Order · Obtain Pricing ·
 - · Check Order Status ·
- (800) 483-5523
- (608) 277-2066 facsimile
- tekladinfo@harlan.com



Helping you do research better