# TD.94059

## Purina #5015, Cocoa Butter, Cholesterol, etc.

Formula		g/Kg
5015, PMI Mouse Diet		750.0
Casein		75.0
Dextrose, monohydrate		30.0
Sucrose		16.25
Dextrin		16.25
Cocoa Butter		75.0
Cholesterol		12.5
Cellulose		12.5
Mineral Mix, AIN-76 (1709	915)	8.75
Vitamin Mix, Teklad (4006	60)	2.5
Choline Chloride		1.25
	051/TD.90221 (Paigen diet) wi	
Colocial Nutriant Int	formet:on1	
Selected Nutrient Inf	% by weight	% kcal from
Protein	19.7	20.4
Carbohydrate <sup>2</sup>	41.2	42.7
Fat	15.8	36.9
Kcal/g 3.9		
	ingredient analysis or manufacturer da lydrate	ata
Teklad Diets are designed & manufactured for research purposes only.		
Speak With A Nutritionist		
• (800) 483-5523		
• askanutritionist@harlan.com		

# Key Features

- Mix of Natural & Refined Ingredients
- Atherogenic
- Cholesterol
- Excludes Cholate

## **Key Planning Information**

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Lead time:
  - · 2 weeks non-irradiated
  - · 4 weeks irradiated

## **Product Specific Information**

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

#### **Options** (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

### **International Inquiry**

- ·Outside U.S.A or Canada·
- askanutritionist@harlan.com

### Place Your Order (U.S.A & Canada)

- Place Order · Obtain Pricing · · Check Order Status ·
- (800) 483·5523
- (608) 277-2066 facsimile

harlan™

• tekladinfo@harlan.com

Helping you do research better

www.harlan.com

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220