# CA.40077 Vitamin Mix AIN-76A

Formula	g/Kg
Thiamin HCI	0.6
Riboflavin	0.6
Pyridoxine HCI	0.7
Niacin	3.0
Calcium Pantothenate	1.6
Folic Acid	0.2
Biotin	0.02
Vitamin B <sub>12</sub> (0.1% in mannitol)	1.0
Vitamin A Palmitate (500,000 IU/g)	0.8
Vitamin E, DL-alpha tocopheryl acetate (500 IU/g)	10.0
Vitamin D <sub>3</sub> , cholecalciferol (400,000 IU/g in sucrose)	0.25
Vitamin K, MSB complex	0.15
Sucrose, fine ground	981.08

### Footnote

Designed for use at 1% diet level (10 g/Kg). This vitamin mix was designed without a choline source because choline bitartrate is listed as a separate item in the formula of the AIN-76A purified diet. References: Second Report of ad hoc Committee on Standards for Nutritional Studies (1980). J. Nutrition 110, 1726. Report of the American Institute on Nutrition ad hoc Committee on Standards for Nutritional Studies (1977). J. Nutrition 107, 1340-1348.

Teklad Mixes are designed & manufactured for research purposes only.

### Speak With A Nutritionist

• (800) 483-5523

www.harlan.com

• askanutritionist@harlan.com



## **Key Features**

- Vitamin Mix
- AIN-76A
- Commonly Used
- Rodent

### **Key Planning Information**

- Products are made fresh to order
- Lead time 2 weeks
- Store product at 4°C or lower
- Use within 6 months (applicable to most mixes)
- Box labeled with product name, manufacturing date, and lot number

### **Product Specific Information**

- Powder
- Minimum order 500 gms

#### **Options** (Fees Will Apply)

- Rush order (pending availability)
- Vacuum packaging (1 and 2 Kg)

### **International Inquiry**

- ·Outside U.S.A or Canada ·
- askanutritionist@harlan.com

#### Place Your Order (U.S.A & Canada)

- Place Order · Obtain Pricing · · Check Order Status ·
- (800) 483-5523
- (608) 277-2066 facsimile
- tekladinfo@harlan.com

Helping you do research better

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220