TD.95027 Low Calcium Diet (0.4% P)

Formula	g/Kg
Casein	200.0
L-Cystine	3.0
Sucrose	342.188
Corn Starch	320.0
Soybean Oil	60.0
Cellulose	40.0
Mineral Mix, Ca-P Deficient (79055)	13.37
Potassium Phosphate, monobasic	11.43
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.012

Key Features

- Purified Diet
- Calcium
- Phosphorus

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
 More frequent replacement may be advised
- Lead time:
 - · 2 weeks non-irradiated
 - · 4 weeks irradiated



Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Footnote Options (Fees Will Apply)

Casein contributes approximately 0.01% calcium or less. Diet contains about 0.4% phosphorus and 2200 IU vitamin D/kg diet. TD.97191 is a possible ingredient matched control with ~0.6% calcium. For questions about this or other calcium adjusted diets askanutritionist@harlan.com.

- Rush order (pending availability)Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Selected	Nutrient I	Information ¹

	% by weight	% kcal from
Protein	17.7	18.6
Carbohydrate	63.4	66.7
Fat	6.2	14.7


Kcal/q 3.8

Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523 option 4
- askanutritionist@harlan.com







www.harlan.com

- · Obtain Pricing · Check Order Status ·
- tekladinfo@harlan.com

Contact Us

• (800) 483-5523 option 3

International Inquiry (Outside USA or Canada)

askanutritionist@harlan.com

Place Your Order (USA & Canada)

- · Please Choose One ·
- https://webapps.harlan.com/teklad-orders/
- tekladinfo@harlan.com
- (800) 483-5523 option 3
- (608) 277-2066 facsimile

Helping you do research better