**TD.95125** Iodine Defic. Diet (0.15% PTU)

Formula	g/Kg
Corn, low iodine	773.0
Wheat Gluten, low iodine	194.0
L-Lysine HCI	1.5
Calcium Carbonate	10.0
Sodium Chloride	10.0
Vitamin Mix, Teklad (40060)	10.0
Propylthiouracil (PTU)	1.5

# **Key Features**

- Mix of Natural & Refined Ingredients
- Modified Remington Diet
- Propylthiouracil (PTU)
- Rodent

# **Key Planning Information**

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
   More frequent replacement may be advised
- Lead time:
  - · 2 weeks non-irradiated
  - · 4 weeks irradiated



# **Product Specific Information**

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

#### **Footnote**

This formula is a modification of TD.95007 to add 0.15% PTU. Each lot of corn and wheat gluten is tested to minimize levels of background iodine. TD.97350 is a possible control without PTU with iodine added at 0.7 mg/Kg from potassium iodate. For questions about this or other iodine adjusted diets askanutritionist@harlan.com.

## Selected Nutrient Information<sup>1</sup>

	% by weight	% kcal from
Protein	20.1	25.4
Carbohydrate	51.4	64.8
Fat	3.5	9.8


Kcal/q 3.2

<sup>1</sup>Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

### **Speak With A Nutritionist**

- (800) 483-5523 option 4
- askanutritionist@harlan.com







www.harlan.com

#### Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

#### **Contact Us**

- · Obtain Pricing · Check Order Status ·
- tekladinfo@harlan.com
- (800) 483-5523 option 3

### International Inquiry (Outside USA or Canada)

askanutritionist@harlan.com

### Place Your Order (USA & Canada)

- · Please Choose One ·
- https://webapps.harlan.com/teklad-orders/
- tekladinfo@harlan.com
- (800) 483-5523 option 3
- (608) 277-2066 facsimile

Helping you do research better