TD.85419

Zinc Deficient Diet

Formula	g/Kg
Egg White Solids, spray-dried	200.0
Dextrose, monohydrate	634.2658
Corn Oil	100.0
Cellulose	30.0
Vitamin Mix, Teklad (40060)	10.0
Biotin	0.004
Mineral Mix, Zn Deficient (81264)	25.6902
Chromium Potassium Sulfate, dodecahydrate	0.02
Ethoxyquin, antioxidant	0.02

Key Features

- Purified Diet
- Zinc
- Egg White Solids
- Rodent

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
 More frequent replacement may be advised
- Lead time:
 - · 2 weeks non-irradiated
 - · 4 weeks irradiated

Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

Contact Us

- · Obtain Pricing · Check Order Status ·
- tekladinfo@harlan.com
- (800) 483-5523 option 3

International Inquiry (Outside USA or Canada)

• askanutritionist@harlan.com

Place Your Order (USA & Canada)

- · Please Choose One ·
- https://webapps.harlan.com/teklad-orders/
- tekladinfo@harlan.com
- (800) 483-5523 option 3
- (608) 277-2066 facsimile
- Helping you do research better

Footnote

Egg white solids are used instead of casein to reduce background zinc. The zinc content of the diet ranges from approximately 0.5 - 1.5 ppm. TD.85420 is a possible ingredient matched control with ~50ppm zinc from zinc carbonate. For questions about this or other zinc adjusted diets askanutritionist@harlan.com.

Selected Nutrient Information ¹		
	% by weight	% kcal from
Protein	16.1	16.6
Carbohydrate	58.3	60.2
Fat	10.0	23.2
Kcal/g 3.9		

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523 option 4
- askanutritionist@harlan.com

Harlan Laboratories, Inc. PO Box 44220 Madison, WI 53744-4220





www.harlan.com

Harlan, Harlan Laboratories, Helping you do research better, and the Harlan logo are trademarks and trade names of Harlan Laboratories, Inc

© 2008 Harlan Laboratories, Inc.