TD.86489

# **Diet with Adjusted** Sucrose/Corn Starch

Formula	g/Kg
Casein	207.0
DL-Methionine	3.0
Sucrose	320.99
Corn Starch	320.0
Corn Oil	50.0
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.01

# **Key Features**

- Purified Diet
- Adjusted Carbohydrate
- Sucrose & Corn Starch
- Rodent

# **Key Planning Information**

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week More frequent replacement may be advised
- Lead time:
  - · 2 weeks non-irradiated
  - 4 weeks irradiated



# **Product Specific Information**

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

### **Footnote Options** (Fees Will Apply)

A diet with approximately equal amounts (wt/wt) of sucrose and corn starch. About 53% of the carbohydrate kcal come from sucrose, and 47% from corn starch. The pellets are relatively short due to the 1:1 ratio of corn starch and sucrose.

•	Rush order (pending availability)	
•	rradiation (see Product Specific Information)	1)

Vacuum packaging (1 and 2 Kg)

Selected Nutrient Information <sup>1</sup>				
	% by weight	% kcal from		
Protein	18.3	20.0		
Carbohydrate	61.7	67.3		
Fat	5.2	12.8		
Kcal/g	3.7			
<sup>1</sup> Values are calculated from ingredient analysis or manufacturer data				

Teklad Diets are designed & manufactured for research purposes only.

## Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220 www.harlan.com

### International Inquiry

·Outside U.S.A. or Canada ·

askanutritionist@harlan.com

### Place Your Order (U.S.A. & Canada)

- · Place Order · Obtain Pricing ·
  - · Check Order Status ·
- (800) 483-5523
- (608) 277-2066 facsimile
- tekladinfo@harlan.com

Helping you do research better

