TD.02028

++++ ENVIGO

Atherogenic Rodent Diet

Formula	g/Kg
Casein	195.0
DL-Methionine	3.0
Sucrose	325.46
Corn Starch	150.0
Anhydrous Milkfat	210.0
Cholesterol	12.5
Cholic Acid	5.0
Cellulose	50.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Carbonate	4.0
Vitamin Mix, Teklad (40060)	10.0
Ethoxyquin, antioxidant	0.04

Key Features

- + Purified Diet
- + Modification of TD.88137
- + Cholesterol
- + Cholic Acid

Key Planning Information

- + Products are made fresh to order
- + Store product at 4°C or lower
- + Use within 6 months (applicable to most diets)
- + Box labeled with product name, manufacturing date, and lot number
- + Replace diet at minimum once per week

 More frequent replacement may be advised
- + Lead time:
 - · 2 weeks non-irradiated
 - · 4 weeks irradiated

Product Specific Information

- + 1/2" Pellet or Powder (crumbly)
- + Minimum order 3 Kg
- + Irradiation available upon request

Options (fees will apply)

- + Rush order (pending availability)
- + Irradiation (see Product Specific Information)
- + Vacuum packaging (1 and 2 Kg)

This is a modification of diet TD.88137 to increase added cholesterol to 1.25% (total in the diet is 1.3%) and add 0.5% cholic acid. TD.96121 is a version without cholic acid. For questions about this diet or modifications askanutritionist@envigo.com

Selected Nutrient Information¹

% by weight	% kcal from
17.3	15.4
46.9	41.9
21.2	42.6
	17.3 46.9

Kcal/g 4.5

Footnote

¹ Values are calculated from ingredient analysis or manufacturer data

Speak With A Nutritionist

- + (800) 483-5523
- + askanutritionist@envigo.com

Teklad diets are designed & manufactured for research purposes only.

Contact Us

Obtain pricing · Check order status

- + teklad@envigo.com
- + (800) 483-5523



International Inquiry (outside USA or Canada)

+ askanutritionist@envigo.com

Place Your Order (USA & Canada)

Please Choose One

- + envigo.com/teklad-orders
- + tekladorders@envigo.com
- + (800) 483-5523
- + (608) 277-2066 facsimile

09/22/15