TD.06415 Adjusted Calories Diet (45/Fat)

Formula	g/Kg
Casein	245.0
L-Cystine	3.5
Corn Starch	85.0
Maltodextrin	115.0
Sucrose	200.0
Lard	195.0
Soybean Oil	30.0
Cellulose	58.0
Mineral Mix, AIN-93G-MX (94046)	43.0
Calcium Phosphate, dibasic	3.4
Vitamin Mix, AIN-93-VX (94047)	19.0
Choline Bitartrate	3.0
Red Food Color	0.1

Footnote

Approx. 45% of total calories come from fat. Designed with similarities to Research Diets, Inc. formula D12451. For the series TD 06414-TD 06416. Approximate fatty acid profile (% of total fat): 36% saturated, 47% monounsaturated, 17% polyunsaturated.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	21.7	19.0
Carbohydrate	41.4	36.2
Fat	22.8	44.8
Kcal/g 4.6		

¹ Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220 www.harlan.com

Key Features

- Purified Diet
- Diet Induced Obesity
- High Fat

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week More frequent replacement may be advised
- Lead time:
 - · 2 weeks non-irradiated
 - 4 weeks irradiated

Product Specific Information

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- Outside U.S.A. or Canada
- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

 Place Order · Obtain Pricing · · Check Order Status ·

(800) 483-5523

harlan™

- (608) 277-2066 facsimile
- tekladinfo@harlan.com
- Helping you do research better