## **TD.89247** 60% Fructose Diet

Formula	g/Kg
Casein	207.0
DL-Methionine	3.0
Fructose	600.0
Lard	50.0
Cellulose	79.81
Mineral Mix, Rogers-Harper (170760)	50.0
Zinc Carbonate	0.04
Vitamin Mix, Teklad (40060)	10.0
Green Food Color	0.15

# **Key Features**

- Purified Diet
- Fructose
- Color Coded
- Rodent

## **Key Planning Information**

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week More frequent replacement may be advised
- Lead time:
  - · 2 weeks non-irradiated
  - 4 weeks irradiated



## **Product Specific Information**

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

#### **Footnote**

This formula is a modification of diet TD 78463 to increase zinc levels. The vitamin mix contributes 0.4% carbohydrate from the corn starch diluent. Fructose is hygroscopic, and the diet will become tacky after a few days of exposure to the ambient environment.

## **Options** (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

## Selected Nutrient Information<sup>1</sup> % by weight

% kcal from Protein 18.3 20.2 Carbohydrate 60.4 66.8 Fat 5.2 12.9

Kcal/g 3.6

<sup>1</sup>Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

#### Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com



Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220 www.harlan.com

#### **International Inquiry**

Outside U.S.A. or Canada

askanutritionist@harlan.com

#### Place Your Order (U.S.A. & Canada)

- · Place Order · Obtain Pricing ·
  - · Check Order Status ·
- (800) 483-5523
- (608) 277-2066 facsimile
- tekladinfo@harlan.com

Helping you do research better

