TD.95217 Adjusted Fat Diet

Formula	g/Kg
Casein	227.7
DL-Methionine	3.3
Sucrose	150.0
Corn Starch	205.24
Maltodextrin	100.0
Vegetable Shortening, hydrogenated (Primex)	106.3
Anhydrous Milkfat	40.0
Soybean Oil	40.0
Cellulose	72.5
Mineral Mix, AIN-93G-MX (94046)	38.5
Calcium Phosphate, dibasic	2.42
Vitamin Mix, Teklad (40060)	11.0
Choline Bitartrate	3.0
Ethoxyquin, antioxidant	0.04

Footnote

A diet used to induce obesity. Approximate fatty acid composition (% of total fat): 34% saturated, 18% trans, 29% monounsaturated (cis), 19% polyunsaturated (cis). Cited in Exp Clin Endocrinol Diabetes 2001; 109: 307-319.

Selected Nutrient Information¹

	% by weight	% kcal from
Protein	20.1	18.8
Carbohydrate	44.3	41.4
Fat	18.9	39.8
Kcal/g 4.	3	

¹Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

Speak With A Nutritionist

- (800) 483-5523
- askanutritionist@harlan.com

Harlan Laboratories · PO Box 44220 · Madison, WI 53744-4220 www.harlan.com

Key Features

- Purified Diet
- Diet Induced Obesity
- Trans Fatty Acids

Key Planning Information

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week More frequent replacement may be advised

• Lead time:

- · 2 weeks non-irradiated
- 4 weeks irradiated

Product Specific Information

- 1/2" Pellet or Powder (crumbly)
- Minimum order 3 Kg
- Irradiation available upon request

Options (Fees Will Apply)

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

International Inquiry

- ·Outside U.S.A. or Canada·
- askanutritionist@harlan.com

Place Your Order (U.S.A. & Canada)

 Place Order · Obtain Pricing · · Check Order Status ·

(800) 483-5523

harlan™

- (608) 277-2066 facsimile
- tekladinfo@harlan.com

Helping you do research better

Harlan, Harlan Laboratories, Helping you do research better, and the Harlan logo are trademarks and trade names of Harlan Laboratories, Inc. © 2008 Harlan Laboratories. Inc.