**TD.96348** Diet (20% Lactose, 2.0% Ca, 1.25% P)

Formula	g/Kg
Casein, "Vitamin-Free" Test	192.0
DL-Methionine	3.0
Sucrose	264.78
Lactose, monohydrate	200.0
Corn Starch	150.0
Corn Oil	50.0
Cellulose	50.0
Vitamin Mix, Teklad (40060)	10.0
Mineral Mix, AIN-76 (170915)	35.0
Calcium Phosphate, dibasic	31.17
Calcium Carbonate	14.05

# Key Features

- Purified Diet
- Lactose
- Calcium
- Phosphorus

# **Key Planning Information**

- Products are made fresh to order
- Store product at 4°C or lower
- Use within 6 months (applicable to most diets)
- Box labeled with product name, manufacturing date, and lot number
- Replace diet at minimum once per week
   More frequent replacement may be advised
- Lead time:
  - · 2 weeks non-irradiated
  - · 4 weeks irradiated



## **Product Specific Information**

- 1/2" Pellet or Powder (free flowing)
- Minimum order 3 Kg
- Irradiation available upon request

## Footnote Options (Fees Will Apply)

Referred to in the literature as a "rescue diet" for mice lacking a functional vitamin D receptor. Contains approximately 2200 IU vitamin D/kg diet, along with 20% lactose, 2% calcium, and 1.25% phosphorus. TD.94112 is a version with green color coding. TD.87095 is a version without vitamin D that is colored brown. For questions about this or other vitamin and mineral adjusted diets askanutritionist@harlan.com.

adjusted diets askanutritionist@harlan.co
Selected Nutrient Information<sup>1</sup>

	% by weight	% kcal from
Protein	17.6	19.8
Carbohydrate	59.8	67.5
Fat	5.0	12.7

Kcal/a 3.5

<sup>1</sup>Values are calculated from ingredient analysis or manufacturer data

Teklad Diets are designed & manufactured for research purposes only.

#### **Speak With A Nutritionist**

- (800) 483-5523 option 4
- askanutritionist@harlan.com

Harlan Laboratories, Inc. PO Box 44220 Madison, WI 53744-4220





www.harlan.com

- Rush order (pending availability)
- Irradiation (see Product Specific Information)
- Vacuum packaging (1 and 2 Kg)

#### **Contact Us**

- Obtain Pricing · Check Order Status ·
- tekladinfo@harlan.com
- (800) 483-5523 option 3

#### International Inquiry (Outside USA or Canada)

askanutritionist@harlan.com

#### Place Your Order (USA & Canada)

- · Please Choose One ·
- https://webapps.harlan.com/teklad-orders/
- tekladinfo@harlan.com
- (800) 483-5523 option 3
- (608) 277-2066 facsimile

Helping you do research better